

# 'Y' Cable Adaptor



**PRODUCT CODE: MED-NC-Y-CABLE-T(X)**

Indigo Care provides a special "Y" Cable Adaptor to be used with all of our Fall Prevention sensor mats and pads. This comes as a standard item when you purchase a Floor Mat or Floor Cushion Mat. The "Y" Cable is sold as a separate item with our chair, bed and wet mats. The "Y" Cable will allow you to connect the sensor mats and pads to your existing nurse call pendant system. The resident will then have a pendant and a mat/pad connected to the one nurse call point. There are three types of "Y" cables to choose from. This choice will depend on your existing nurse call system. The cord is 3 meters long which will allow for greater coverage as well.

## **TYPE 1: MED-NC-Y-LEAD-T1**

- **TYPE 1 : Stereo Plug**
- **Suitable for 1 and 2 button nurse pendant with 6.3mm stereo**

## **TYPE 2: MED-NC-Y-LEAD-T2**

- **TYPE 2: Mono Plug.**
- **Suitable for 1 button nurse pendant with 6.3mm mono plug**

## **TYPE 3: MED-NC-Y-LEAD-T3**

- **TYPE 3: Stereo Plug**
- **Suitable for 1 and 2 button nurse pendant with 6.3mm stereo plug sleeve common**

**PORTABLE**

**EASY TO USE**

**3 TYPES OF "Y" CABLE ADAPTORS, T1, T2, T3**

**USED WITH:  
FLOOR SENSOR MATS  
BED SENSOR MATS  
BED SENSOR PADS  
CHAIR SENSOR MATS  
WET SENSOR MATS  
WET SENSOR PADS**

**INTEGRATES TO MOST NURSE CALL SYSTEMS**

**PLUG & PLAY OPERATION**

**3 METERS IN LENGTH**

**MOULDED BOTH ENDS**

**WARRANTY:  
12 MONTHS**



**HORIZON ELECTRONIC COMPONENTS Pty Ltd**  
2/71 Rushdale Street, Knoxfield, VIC, 3180  
PH +61 3 9763-2988 FAX +61 3 9763-3811  
[www.horizone.com.au](http://www.horizone.com.au)  
ABN 36007302801



## **ORDER INFORMATION**

**Fax or email your purchase orders to  
03-9763-3811  
Or email [sales@horizone.com.au](mailto:sales@horizone.com.au)**