MULTI PURPOSE ALARM MANAGEMENT:

RECORDABLE VOICE
MESSAGE FOR ALARM

NORMALLY OPEN SWITCH

NORMALLY CLOSE SWITCH

**MAGNETIC SWITCH** 

REDUCE THE RISK OF FALLS & INJURIES

**PORTABLE** 

**EASY TO USE** 

EASY TO MOVE FROM ROOM TO ROOM

INTEGRATES TO MOST NURSE CALL SYSTEMS

PLUG & PLAY OPERATION

**DURABLE & RELIABLE** 



## **Voice Alert Sensor Monitor**



PRODUCT CODE: MED-NC-TS0893

Indigo Care provides the Voice Alert Monitor that is designed to be used in conjunction with our floor mats and bed mats/pads. This is a multi-function fall management alarm which can be used as a monitor with magnetic switch, floor sensor mat and bed/chair sensor mat. There is the choice of voice or four different alarm tones which will sound when a resident leaves their bed/chair or steps on the floor mat or when the magnetic switch is detached (when magnetic switch is in use). Record your own alert message or instructions in a voice familiar to the resident that will play in conjunction with the alarm. The monitor can connect and alarm through your current nurse call system. There are two alarm sounds & a Low or High Volume setting. Alternatively, this alarm sound can be put on silent mode and used with your nurse call system.

## **KEY SPECIFICATIONS:**

- Versatile mounting base
- Choice of voice or 4 different alarm tones
- Record your own alert message or instructions in a voice that is familiar to the resident.
- Reset button
- 2 alarm sounds High & Low
- Silent mode (with flashing light) for use with remote alarm signaling Nursecall systems
- Easy Clean Pull Cord
- Handy compartment to store pull cord at back of monitor
- Power: Battery Operated 4 x AA (not included)
- Low Battery warning
- Warranty: 12 months

