#### **PORTABLE**

**EASY TO USE** 

3 TYPES OF "Y" CABLE ADAPTORS, T1, T2, T3

USED WITH:
FLOOR SENSOR MATS
FLOOR CUSHION MATS
BED SENSOR MATS
BED SENSOR PADS
CHAIR SENSOR MATS
WET SENSOR MATS
WET SENSOR PADS

INTEGRATES TO MOST NURSE CALL SYSTEMS

PLUG & PLAY OPERATION

3 METERS IN LENGTH

**MOULDED BOTH ENDS** 

WARRANTY: 12 MONTHS



## 'Y' Cable Adaptor



# PRODUCT CODE: MED-Y-LEAD-T(X) (TGA Registration: 223902)

Indigo Care provides a special "Y" Cable Adaptor to be used with all of our Fall Prevention sensor mats and pads. This comes as a standard item when you purchase a Floor Mat or Floor Cushion Mat. The "Y" Cable is sold as a separate item with our chair, bed and wet mats. The "Y" Cable will allow you to connect the sensor mats and pads to your existing nurse call pendant system. The resident will then have a pendant and a mat/pad connected to the one nurse call point. There are three types of "Y" cables to choose from. This choice will depend on your existing nurse call system. The cord is 3 meters long which will allow for greater coverage as well.

#### **TYPE 1: MED-NC-Y-LEAD-T1**

- TYPE 1: Stereo Plug
- Suitable for 1 and 2 button nurse pendant with 6.3mm stereo plug tip common

#### TYPE 2: MED-NC-Y-LEAD-T2

- TYPE 2: Mono Plug
- Suitable for 1 button nurse pendant with 6.3mm mono plug

### **TYPE 3: MED-NC-Y-LEAD-T3**

- TYPE 3: Stereo Plug
- Suitable for 1 and 2 button nurse pendant with 6.3mm stereo plug sleeve common

Horizon Electronic Components T/as INDIGO CARE 2/71 Rushdale Street, Knoxfield, VIC, 3180 PH +61 3 9763-2988 FAX +61 3 9763-3811 www.indigocare.com.au ABN 36007302801

ORDER INFORMATION
Fax or email your purchase orders to
03-9763-3811
Or email sales@indigocare.com.au