

SUITABLE FOR UNDER
THE MATTRESS

REDUCE THE RISK OF
FALLS & INJURIES

LOW PROFILE,
NON-SLIP MATERIAL

PORTABLE

EASY TO USE

EASY TO MOVE FROM
ROOM TO ROOM

SOLVENT WIPE-
DISINFECTANT

INTEGRATES TO
MOST NURSE CALL
SYSTEMS

PLUG & PLAY
OPERATION

DURABLE & RELIABLE

COST EFFECTIVE

LONG PRESSURE MAT
(760mm X 110mm)



Under Mattress Sensor Pad



PRODUCT CODE: MED-B365-UM-4
(TGA Registration: 223902)

Indigo Care provides a pressure sensitive Under Mattress Sensor Mat that is designed for UNDER MATTRESS to immediately alert staff should a resident leaves the bed. This is most useful to residents who are at risk of falling out of their bed, or those who require assistance in arising from the bed. The sensor mat may not be suitable for some mattresses eg. air mattress. The bed sensor mat is connected to your existing nurse call system via the monitor and a Y-Cable. Once the resident removes themselves from the mat and tries to leave the bed, the sensor mat will activate the monitor's alarm and alert the staff. There are a set of instructions printed on the actual mat. Eyelets at the end of the mat are for easy storage.

Manufacturer/ Supplier recommendation:

Check the sensor mat is compatible for the mattress before purchased. Get a trial sample from Indigo Care

KEY SPECIFICATIONS:

- **Dimensions: 760mm X 110mm**
- **Connection: RJ11 on 500mm curly cord**
- **Electrical Resistance <100 Ohm**
- **White PVC, Thickness 0.20mm**
- **Warranty: 12 months**

PRODUCT ACCESSORIES:

- **Y Cable**
- **Anti-temper Monitor (use 4 x AA batteries-not included)**
- **AC Adaptor (optional)**

Horizon Electronic Components T/as INDIGO CARE
2/71 Rushdale Street, Knoxfield, VIC, 3180
PH +61 3 9763-2988 FAX +61 3 9763-3811
www.indigocare.com.au
ABN 36007302801



ORDER INFORMATION

Fax or email your purchase orders to
03-9763-3811
Or email sales@indigocare.com.au