


Wet Sensor Pad



PRODUCT CODE: MED-NC-WP-01

Indigo Care provides a Wet Sensor Pad that is designed to immediately alert staff should a resident wet their bed or chair. This is most helpful to residents with incontinence issues. The pad works in conjunction with the Deluxe Wet Sensor Monitor. The Wet Sensor Pad is placed on the bed, directly under the bed linen. The Deluxe Wet Sensor Monitor is connected to the wet sensor pad and then, with a Y cable, the pad is connected to your existing nurse call system via the monitor. If the resident were to have an "accident", the pad's sensor strips will detect the moisture and the alarm will sound alerting the staff. The pad is waterproof and easy to wipe down and clean. A set of instructions is printed on the back of the pad.

**REDUCE THE RISK OF
INCONTINENCE
ACCIDENTS**

**LOW PROFILE,
NON-SLIP MATERIAL**

PORTABLE

EASY TO USE

**EASY TO MOVE FROM
ROOM TO ROOM**

WATERPROOF

WASHABLE

**INTEGRATES TO MOST
NURSE CALL SYSTEMS**

**PLUG & PLAY
OPERATION**

DURABLE & RELIABLE

COST EFFECTIVE

**LONG WET SENSOR
MAT
(760mm X 110mm)**

KEY SPECIFICATIONS:

- **Dimensions: 760mm X 110mm**
- **Connection: Curly Cord on RJ11 (4P4C)**
- **Resistance <10k Ohm**
- **IP Rating : 66**
- **White & Green PVC, Thickness 0.20mm**
- **Cord Length 2 Meters**
- **Warranty: 12 months**

PRODUCT ACCESSORIES:

- **Y Cable** 
- **Deluxe Wet Sensor Monitor** 
- **AC Adaptor (optional)**

indigo

care

HORIZON ELECTRONIC COMPONENTS Pty Ltd
2/71 Rushdale Street, Knoxfield, VIC, 3180
PH +61 3 9763-2988 FAX +61 3 9763-3811
www.horizone.com.au
ABN 36007302801



ORDER INFORMATION

Fax or email your purchase orders to
03-9763-3811
Or email sales@horizone.com.au