

# Pull Cord Pendant

MED-NC-PCP-V2.0



## Key Specifications:

- Call plate is ceiling mounted
- Mounting Plate – White & Plastic
- Light Indicator – Red LED
- Cord is 2.5 metres & Nylon
- Cord breaking strength – 10 kg
- Call Plate Dimensions:  
110 mm (L) x 67 mm (W) x 14 mm (D)
- Momentary Action Switch
- Protection against moisture & dust
- Optional: Variable Cord Length

Indigo Care has designed the Pull Cord Pendant with the elderly resident in mind. The call plate is mounted on the ceiling so the resident will have easy access to the pull cord to initiate a call. The cord length can be variable to suit different ceiling heights. There is a red LED for call activity. The switch is momentary action and is reset when the cord is released.

Part Number: MED-NC-PCP-V2.0

*Easy to Use*

*Cost Effective*

*Integrates to most  
Nurse Call Systems*



HORIZON ELECTRONIC COMPONENTS Pty Ltd  
2/71 Rushdale Street, Knoxfield, VIC, 3180  
PH +61 3 9763-2988 FAX +61 3 9763-3811  
[www.horizone.com.au](http://www.horizone.com.au)  
ABN 36007302801



## ORDER INFORMATION

Fax or email your purchase orders to  
**03-9763-3811**

Or email [sales@horizone.com.au](mailto:sales@horizone.com.au)